



# Multi Valve Combinations

## Flexibility

**M**ulti combination diverter style diaphragm valves have become a cost effective solution to overcrowding in complex piping layouts. Our customers want the same quality diverter valve they have come to expect from Aquasyn's other products.

**W**e can provide 3, 4, and 5 way diverters, for both horizontal and vertical placement. Aquasyn's actuated and manual valves offer maximum drainability and process flexibility.



## Features

- Machined block construction
- Compact solutions to overcome crowded multiple valve combinations
- Decreased dead legs to remedy complex drainage problems
- Less system welds
- Optimum drainability
- Quick and effective sterilization

## Custom Valve Design

### Exceptional Capabilities

**C**ustomized and hybrid valve assemblies for unusual applications have been Aquasyn's forte since its inception. We don't think anyone does it quite as well. Call us when you need it done right and on time the first time.



## Unique Applications

- Special end connections, configurations and material
- Innovative designs assist in providing application specific process solutions
- Sole source responsibility from concept through production
- CAE design in state of the art 3-D formats in NX-Ideas, SolidWorks & AutoCAD
- Fully detailed drawings for review and approval
- Rapid customer turnaround with concepts, drawings and production
- Exotic materials a specialty
- Expert fabrication, welding, and polishing by skilled craftsmen
- Advanced machining capabilities optimize manufacture of complex assemblies

